

Faecal Sample Examination

Direct Smear Method

- 1) Place one drop of distilled water on clean slide**
- 2) Put a small amount of faecal sample**
- 3) Mix it uniformly using smooth glass rod or match stick**
- 4) Place a coverslip and examine under low power**

(It fails to detect low grade infections)

Flotation Method

- 1) Mix thoroughly 2gm of faeces with 20ml of saturated salt (264gm/lit) or magnesium sulphate (410gm/lit) solution in beaker**
- 2) Pour the suspension in the centrifuge tube upto brim through strainer or muslin cloth**
- 3) Put coverslip touching the surface of fluid**
- 4) Allow it to stand for 20-30 minutes**
- 5) Take off coverslip, put it on slide and examine under low power**

(It is suitable for nematode eggs and oocysts)

Sedimentation Method

- 1) Mix thoroughly 2gm of faeces with 20ml of distilled water in beaker**
- 2) Pour the suspension in the centrifuge tube upto brim through strainer**
- 3) Centrifuge tube at 1000 rpm for 2-3 minutes**
- 4) Discard the supernatant and take 1-2 drops of sediment on slide**
- 5) Put coverslip and examine under low power**

(It is suitable for trematode/fluke eggs)