

Haemoglobin Estimation

[Sahli's hemoglobinometer method]

- 1. N/10 HCl is put in the hemoglobinometer upto the mark 2**
- 2. Draw fresh blood upto 20 cumm mark in the pipette**
- 3. Immediately mix it with N/10 HCl**
- 4. Allow it to stand in the comparator for 10 min. for the conversion to acid hematin**
- 5. Add distilled water after 30 min. drop by drop and mix with stirring rod till the colour matches well with fixed colour in the comparator**
- 6. The mark tallying with the upper level of diluted acid hematin denotes level of Hb**
- 7. It is expressed in terms of g/100 ml**