

# **Toxicological Techniques**

## **For Nitrate/Nitrite Poisoning**

### **Diphenylamine Blue Test:-**

- 1) Place 1-2 drops of plasma/urine/rumen fluid on glass slide**
- 2) Put 1-2 drops of DPB reagent by the side of same and allow to mix**
- 3) Development of Blue colour indicates positive test**

## **For HCN Poisoning**

### **Picrate Paper Test:-**

- 1) Put few drops of Rumen fluid on Picrate paper**
- 2) Change of paper colour from Yellow to Brown/Red indicates presence of Cyanide**